

Measurement of BMDBM(Butyl Methoxydibenzoylmethane), OMC(Ethylhexyl methoxycinnamate)

[SJ-MM-026]

2014-07-16

REV.2

Quality Assurance

Test procedure

1. Put the sample to 10ml volume flask.
2. Add 10ml of methanol and sonicate the flask for 1hr.
3. Centrifuge (3,500rpm, 10min) & Syringe filter (0.45um) to remove impurities.
4. Analysis
5. Analysis condition

Apparatus Model	Agilent 1200 HPLC
Column	ZORBAX Eclipse XDB-C18 250 × 4.6mm I.D. 5um ZORBAX Eclipse XDB-C18 12.5 × 4.6mm I.D. 5um
Mobile phase	0.1% TFA in water : Acetonitrile = 10 : 90
Detector	UV 358 nm(BMDBM), UV 310 nm(OMC),
Column oven temp	35°C
Flow rate	0.5ml/min
Injection volume	1ul



HPLC : High Performance Liquid Chromatography