

Measurement of Particle Size

[SJ-MM-002]

2010-03-30

Supersedes prior issues

Quality Assurance

Test procedure

1. Clean Particle Size Analyzer(① Coulter, LS13320 ② Malvern Mastersizer 2000) vessel.
2. Put a sample(powder, liquid) into 200ml IPA Solution, and stir well.
3. Sonicate the sample about 2 ~ 3 minutes.
4. Measure Particle Size as described in Manual.

Equipment

① Coulter, LS13320



② Malvern Mastersizer 2000

