

<b>Measurement of SPF</b>	<b>[SJ-MM-006]</b>
<b>In Vitro Sun Protection Factor(SPF) Testing Method</b>	2005-05-25 Supersedes prior issues
	<b>Quality Assurance</b>

<p><b><u>Test procedure</u></b></p> <ol style="list-style-type: none"><li>1. Tape substrate : Traspore™ tape, 3M</li><li>2. Sample preparation : Vaseline 2g + powder 1g</li><li>3. Sample amount : 1mg of powder/cm<sup>2</sup> ⇒spread 3mg of mixture(Vaseline 2g + powder 1g) on 1 cm<sup>2</sup> of Traspore™ tape</li><li>4. Measure 12 points of the tape, drop out highest and lowest abnormal data and have average of the other left data.</li></ol>
---

<p><b><u>Equipment</u></b></p> <p>Optometrics 290</p>
---